GARANCIJA

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WARRANTY

- The product will work faultlessly if used in accordance with operating instructions.
- Service for maintenance is provided by the manufacturer.
- The manufacturer is committed to repair any deficiencies and breakdowns of the product - resulting from its normal application - within 45 days from the date the product has entered the repair shop. If the product cannot be repaired the manufacturer is committed to replace the faulty product with a faultless one.
- This warranty repair does not include component parts which are subject to normal wear, e.g. packing elements, bactericidal filters and other filters, coupling pipes and mechanical damage.
- This warranty is valid if the warranty card has been certified by the date of purchase and signed and stamped by the salesman.
- When bringing forward a claim please include the product and the certified warranty card.
- The manufacturer will cover the expenses contracted by transportation of faulty product within warranty period to the amount of valid bus, railroad or postal rates upon submission of receipts.
- If the repair shop of the manufacturer does not fix minor faults the warranty will be extended. In case the product has been replaced or only a part of the product has been repaired the warranty starts for that particular part only.
- This warranty becomes void if repair work has been done by persons not authorized by the manufacturer.
- This warranty is valid for the period of one year from the date of purchase

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GB
Dear customer!
Thank you for
purchasing our
product. We would,
however, like to ask
you to read the
operating instructions
very carefully prior to
application of the
product

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DISINFECTION IN LIQUID DISINFECTANTS

Only humidifier component
parts can be sterilized in an
autoclave at 120°C to 134°C!
All other parts of the
flowmeter must be cleaned
with dry cloth only! Hose
barb for humidifier ID
1140015 is not intended for
sterilization in autoclave at
temperature 121-134°C.
Hose barb can be sterilized
with plasma method.

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WARNING!

- Device shall be use for dosing of medical gases only!
- Handling and servicing requires special attention to the cleanliness of hands, tools and environment!
- Device can be used by properly trained medical staff only!
- Greasing of any component part is strictly prohibited!
- Due to the device working at the pressure of 5 bar utmost!
- Caution is required!
- Changes in inlet pressure can affect the accuracy of gas flow measurement to a great degree!
- For humidifying of gases only distilled water can be used!
- When the humidifier is filled with the water, do not turn the flowmeter upside down! In such case the water could enter into the measuring tube of the flowmeter.
- The product shall be stored in dry places only.
- When changing connector on the Flowmeter never use connector of different gas as indicated on the flowmeter.

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

GB
MAINTENANCE AND
SERVICING

Disconnect the flowmeter from the gas supply before servicing

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Kontrola na 2 leti
Po vsakem popravilu

- Tehnični pregled

Prescribed control before every use
- Functional checkup

Control every 2 years, after each maintenance

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GB
INSTRUCTIONS FOR
SAFE USE

TECHNICAL DATA
Drive gas: oxygen, compressed air
Inlet pressure (p1): 3,5-5,5 bar
Flow rate: 0-15 l/min
Flowmeter is calibrated at 21°C

DESCRIPTION
Flowmeter is designed to add gas for inhalation or insufflation. Regarding the type of gas treatment there are so called direct variants available, which are with an appropriate connector that is attached to an outlet of medical gases and so-called rail variants, where the flowmeter is attached to a rail system and the gas source is brought to the device through a flexible tube with a connector. A flowmeter without humidifier doses dry gas, but a humidifier which humidifies gas with distilled water can be added. Gas enters through the atomizer pout and sprays from the container; due to the minimal size (1µm) of its particles the aqueous aerosol penetrates deep into respiratory tract where most of the water vaporizes which is why the gas is saturated by water vapour so that mucous membrane cannot dry out. The quantity of dispersed water at a flow rate of 15l/min is 36g/h. For humidifying of gas only the distilled water can be used. Regarding the quantity of gas treatment, it can be used single or double rail mounted flowmeters. The advantage of

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
double flowmeter is that two patients can be simultaneously treated from one outlet.

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SL Tirnični priključek z privijalnim gumbom-aluminij


GB Rail brackets with knob – aluminium

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INHALATION

Through a hole in the mask a patient can additionally inhale the air from the environment if the adjusted amount of gas is smaller than the minute volume. In this way a step regulation of gas concentration in the inhaled air is possible.

INSUFFLATION

Catheter is introduced through the nose and all the way to the soft palate. In case the gas amount is smaller than the minute volume, the patient inhales the air from the environment through his nostrils whereby the setting of the gas zone is made possible. When the gas flowmeter is not required any longer it is closed by turning the dosing screw to the right.

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GB

Direct flowmeters

The flowmeter is supplied with gas in that the connector is connected to the appropriate medical gas clutch

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GB



Rail mounted flowmeters

Flowmeter should be placed on the wall rail by the bracket. The flowmeter is supplied with medical gas in that the connector with tube is connected into the appropriate medical gas clutch

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