
RVTM3

- F Notice d'instructions
- GB Instructions for use
- D Bedienungsanleitung
- E Folleto de instrucciones
- I Avvertenze e istruzioni
- P Manual de Instruções
- BG Инструкции за експлоатация
- GR Οδηγίες χρήσης
- H Használati útmutató
- LT Naudojimo instrukcija
- MT Ara n-nota ta' struzzjonijiet
- NL Gebruiksaanwijzing
- PL Instrukcja obsługi
- RO Manual de Instrucțiuni
- SK Návod na použitie
- SL Navodila za uporabo
- SE Bruksanvisning
- CZ Návod k použití
- AR إرشادات الاستخدام

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1. Presentation
2. Versions
3. Technical data
4. Instructions for use
5. Symbols
6. Safety advice
7. Cleaning and disinfection
8. Transport and storage
9. Maintenance
10. Guarantee

1. Presentation

Medical device meant to **measure and adjust the vacuum level during medical or surgical suction**. It connects to a wall vacuum supply. The suction controller should be used with a collection jar and a suction hose.

2. Versions

The RVTM3 Suction Controller is available in the following versions:

- Maximum vacuum levels: 0-250, 0-600, 0-1000 mbar/hPa or 0-400, 0-760 mmHg.
- Inlet fitting: Compliant with NF S 90-116 - BS 5682- DIN 13260 - US Ohmeda Diamond - NORDIC – UNI.
- Direct connection or rail mounted (Hose + crimped hose holder + rail clamp).
- Outlet: nipple or nipple + jar.

3. Technical data

- Compliant with ISO 10079-3: 2009.
- Unit: millibar (mbar), hectoPascal (hPa) and millimetre of mercury (mmHg): 10 mbar = 10hPa = 7.5 mmHg.
- Accuracy of measures: 1.6 % of full scale of the vacuum gauge.
- Vacuum gauge adjustable through 90° (from - 45° to + 45°).
- Free rotation of regulation knob at its end-of-travel position (no possible blocking).
- Safety jar adjustable through 90° (from - 45° to + 45°).
- Antibacterial filter: 99.97% efficiency for 0.3 µm particles.
- Unit serial number: laser engraved on the device body.
- 0-250 mbar version has a safety valve.
- Metal connection nipple, fitted directly into the body of the suction controller.
- Maximum flow rate: 110 l/min (for a pump with a flow rate of 160 l/min).
- Operating temperature: 23°C.

4. Instructions for use

- The suction controller must be used exclusively by people who have been trained.
- Influence of the vacuum supply: caution, check the maximum vacuum level of the vacuum supply. The maximum vacuum level of the suction controller should not exceed the maximum value of the supply.

Suction controller

- Make sure the device is undamaged and that the connector complies with the supply connection.
- Make sure the regulation knob is closed.
- Turn the device OFF (red part of the push button visible).
- Connect the device to the wall vacuum outlet.
- Connect the suction controller's outlet (nipple or safety jar) to the collection jar (Ø 9 mm) using a min. 6.3 diameter tubing.
- To start suction, press the red button to make green pushbutton visible).
- Turn the front regulation knob anti-clockwise and adjust the vacuum level by closing the patient suction system.

→ For a temporary stop, use the ON/OFF pushbutton. Vacuum level comes back to its initial value when suction is back on.

→ For a long-term stop, turn off the regulation knob and the ON/OFF pushbutton (suction system will purge).

If vacuum level has been over-adjusted, shut down the regulation knob and purge the device (open the suction system/initiate a flow) before adjusting the level again.

Safety jar

- **Setting up:** by clipping under suction controller:
Align the mark found on the cover of the safety jar with that found on the skirt attached to the body of the suction controller, then turn the jar ninety degrees clockwise.
- **Removing:** turn the safety jar to the left, align the mark found on the cover of the safety jar with that found on the skirt attached to the body of the suction controller, then pull the jar downwards.

Filter

- **Setting up:** push until click-lock.
- **Removing:** pull while operating a rotational movement.
- Throw the filter away in compliance with standard procedures.

<p>Why an antibacterial filter? Aerosol that may be produced during suction might contaminate suction controllers, fittings, wall outlets, pipelines and vacuum pumps. Moreover, without any flow or filter, bacteria might circulate freely in the patient circuit.</p>

5. Symbols

See instructions for use

GB

6. Safety advice

- Always test prior to use.
- Avoid impacts with the meter (needle is fragile).
- Make sure the needle is on 0 when the device is not connected.
- Make sure the device is always fitted with a filter.

7. Cleaning and disinfection

In normal conditions of use, the RVTM3 does not need any disinfection since the safety jar is protected by the filter upfront. **Change the filter after each patient.** (See below for specific conditions).

Suction controller

Use soapy water for the outside of the device. Rinse thoroughly and dry. Do not immerse.

Caution:

- When using decontaminant products, make sure they are compatible with plastic (see technical data).
- Do not use surface decontaminants.

Safety jar

Since the safety jar is protected by an upfront filter, there is no need to sterilize it: the filter just needs to be replaced after each patient or for a long stay patient, when necessary, after inspection of the filter clogging level. However, in case of accidental liquid overflow or perforated filter, the safety jar must be disinfected or sterilized as follows:

1. Remove the safety jar and disassemble the cover from the jar.
2. Immerse safety jar parts in a pre-disinfecting solution (respecting manufacturers' protocols concerning soaking, rinsing and drying time).
3. Clean, rinse thoroughly, soak for 1 minute in clean water and dry.
4. Sterilize or disinfect, respecting manufacturers' protocol.
5. Reassemble each part and screw the safety jar back onto the suction controller.

RVTM3 safety jar withstands autoclave up to 134°C.

8. Transport and storage

- To prevent damage to gauge needle, devices should travel in a horizontal position.
- Storage between -10 to +40°C in a dry and clean place.
- Preserve the device packaging and instructions for use as long as possible.

9. Maintenance

Check every 1 to 3 years depending on use:

- Inspect display (vacuum gauge).
- Inspect seals.

In case of damage or malfunction, ship the device back to the manufacturer or to a registered repairer.


- Once a year, lubricate the seals of the ON/OFF pushbutton with "High Vacuum Grease" silicon grease
To remove the pushbutton, pull the red area.
- Record every maintenance operation.
- To change the vacuum gauge, remove the 2 screws at the rear of the body.

Spare parts	
- Safety jar	Use exclusively spare parts from <u>TECHNOLOGIE MEDICALE</u>
- "High vacuum Grease"	
- 250 / 600 / 1000 mbar and 400 / 760 mmHg vacuum gauge with protective shell	
- Tube of 10 filters	
- Connector	

10. Guarantee

- 1 year warranty (parts and service), except in case of damage or accidents due to carelessness, misuse, bad supervision or maintenance.
- Contractual lifetime: 10 years.

*As per the European Directive 93/42/CEE dated June 14th 1993,
TECHNOLOGIE MEDICALE must be immediately notified of any incident or risk of incident*

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